

The following shows the effective features of the Upper detector (secondary electron detector: SED), Lower detector (LDD) and the various modes. Select the most appropriate mode for your needs.

WD	Selecting energy range (specimens that do not charge)		Measures for charged specimens			
Upper  1.5-5 mm	SEI		SEI		GB	
	SE	BE	E	EV	GB-High	GB-Low
	Secondary electron (Pure SE) Topographic image	Backscattered electron (Pure BE) Composition image	Energy filter  Capture parts of the secondary electron and eliminate unevenness in the image.  <ul style="list-style-type: none"> <li>• Sensitive to localized electric fields in the specimen.</li> <li>• Irradiation current cannot be increased.</li> <li>• Sensitive to current density.</li> <li>• Some specimens cannot be viewed in Fine View.</li> </ul>		(High BIAS) 1.5kV-2kV  GB mode for capturing information on the extremely shallow surface of the specimen. (inclination not allowed)	(Low BIAS) 100-400V  GB mode for eliminating charge (inclination allowed)  <ul style="list-style-type: none"> <li>* Charge can be suppressed by inclining the specimen.</li> <li>• Reduced diffusion domain</li> <li>• Electron extrusion effect due to the electric field.</li> <li>• Due to combined effects such as capturing the high-energy secondary electrons.</li> </ul>
	Sb	Bs	→ This is a compensatory method.		For extremely-low acceleration and high resolution (0.1kV)  Dynamic imaging of topographic images  <Reference> GB is a type of lens based on electromagnetic field superposition.	
					Measures against charge by secondary electrons	Measures against charge by secondary and backscattered electrons
Lower  4-12 mm	LEI		LEI			
					Acquisition of high-energy secondary electron. (Low-angle backscattered electrons) -> Sensitive to localized electric fields on the specimen. The irradiation current (PCc) can be increased.  Slow Scan is possible. Image quality: Composition image + Topographic image <u>Creates a better topographical image than SEI (Upper detector).</u>  Supports the Bs of the upper detector. <ul style="list-style-type: none"> <li>• Measures against charge up to x100K</li> <li>• Effective during analyses.</li> </ul>	Measures against charge by backscattered electrons